



RECEIVED
AUG 05 2002
TECH CENTER 1600/2900

Page A-1

**APPENDIX A - SPECIFICATION/CLAIM AMENDMENTS
INCLUDING NOTATIONS TO INDICATE CHANGES MADE**

**Serial No.: 09/981,286
Docket No.: 265.00260101**

Amendments to the following are indicated by underlining what has been added and bracketing what has been deleted. Additionally, all amendments have been shaded.

In the Specification

The paragraph beginning at page 19, line 2, has been amended as follows:

Cells that are useful in the methods described herein vary depending on the pathogen or toxin used. The cell chosen for use with a particular pathogen or a particular toxin is one for which the pathogen or toxin is cytotoxic, preferably, cytolytic. Which cells are appropriate for use with a particular pathogen or a particular toxin in known to the art, and can be chosen by a person having skill in the art. Examples of some pathogen/cell combinations include *Bacillus anthracis* and the macrophage cell line J774A.1 (Friedlander et al., *Infect. Immun.*, 61, 245-252 (1993), and Hanna et al., *Mol. Biol., Cell*, 3, 1269-1277 (1992)), Rift Valley fever virus and Vero cells, Venezuelan equine encephalitis and Vero cells, and *Rickettsia* spp. and primary chick embryo cells (Walker and Cain, *Lab Invest.*, 43, 388-396 (1980)) and human endothelial cell culture (Walker et al., [*Fed. Proc.*, 40, 72A] *Federation of American Societies for Experimental Biology*, 40(3):776 (1981)).

The paragraph beginning at page 29, line 18, has been amended as follows:

The adaptein library-containing vector DNA is isolated from *E. coli* bacteria by standard techniques and purified by equilibrium density gradient separation using cesium chloride to form the gradient and following standard techniques. The DNA (adaptein library) is dissolved at 1 mg/ml in 1mM EDTA, 10 mM Tris-HCl, pH 8.0 (TE) and stored at -80°C until required. To produce the MLV encoding the adaptein library, HEK 293 human fibroblasts are grown to 80% confluence on plastic plates in Dulbecco's Modified Eagle Medium (DMEM) supplemented with 10% fetal calf serum (culture medium) and with the bacterial antibiotics streptomycin and penicillin. The cells are simultaneously transfected with vector DNA encoding the adaptein library (25 μg), the MLV structural genes gag-pol (25 μg) and the Vesicular Stomatitis virus glycoprotein (VSV-G; 5 μg). Transfection is conducted according to the methods of Chen and Okayama (*Mol. Cell. Biol.*, 7, [2745-5272] 2745-2752 (1987) and used calcium phosphate to form a precipitate with the DNA that is efficiently taken up by this cell type and the encoded genes expressed. The DNA is mixed with 0.1 mM EDTA, 1 mM Tris-HCl, pH 8.0 in purified water to give a total volume of 1.25 ml. To this 1.25 ml of 50 mM BES, 280 mM NaCl, 1.5 mM Na_2HPO_4 , pH 7.0 is added and mixed. Calcium chloride (2 M, 0.16 ml) is added dropwise while mixing and then incubated for 10 minutes at room temperature. This is then applied to cells with 2.5 ml of the mixture being used per 15 cm diameter plastic dish containing 25 ml of culture medium. The following day the culture medium is removed and cells washed with 7 ml of fresh medium and replaced with 20 ml of culture medium. To generate sufficient complexity in the library this procedure is repeated 50 times to yield 1×10^8 transfected cells expressing each library member.